

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

**Recombinant Streptococcus pyogenes serotype M3
Phosphoribosylformylglycinamide cyclo-ligase(purM), Escherichia Coli
Preis auf Anfrage
CSB-EP320467SMV**

Artikelname	Recombinant Streptococcus pyogenes serotype M3 Phosphoribosylformylglycinamide cyclo-ligase(purM), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP320467SMV
Hersteller Artikelnummer	CSB-EP320467SMV
Alternativnummer	CSB-EP320467SMV-1, CSB-EP320467SMV-100, CSB-EP320467SMV-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	P0DD62
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MSEKNAYAKS GVDVEAGYEV VERIKKHVAR TERAGVMGVL GGFGGMFDLS KTGVPKEPVLV SGTG DGVGTKL MLAIKYDKHD TIGQDCVAMC VNDIIAAGAE PLYFLDYIAT GKNNPVKLEE VVSGVAEGCV QAGVALIGGE TAEMPGMYGE DDYDLAGFAV GVAEKSQIID GSKVKEGDIL LGLASSGIHS NGYSLVRRVF ADYTGKELLP ELEGKQLKDV LLEPTRIYVK AA
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.