

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Escherichia coli HTH-type transcriptional activator rhaS(rhaS),Escherichia Coli Preis auf Anfrage CSB-EP482341ENM

Artikelname	Recombinant Escherichia coli HTH-type transcriptional activator rhaS(rhaS),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP482341ENM
Hersteller Artikelnummer	CSB-EP482341ENM
Alternativnummer	CSB-EP482341ENM-1, CSB-EP482341ENM-100, CSB-EP482341ENM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	E. coli
Tag	The tag type will be determined during production process.
UniProt	B7L9G1
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTVLHSVDFP PSGNASVAIE PRLPQADFPE HHHDFHEIVI VEHGTGIHVF NGQPYTITGG TVCFVRDHR HLYEHTDNLC LTNVLYRSPD RFQFLAGLNQ LLPQELDGQY PSHWRVNHVS LQQVRQLVAQ MEQQEGENDL PSTASREILF MQLLLLLRKS SLQENLENSA SRLNLLLAWL EDHFADEVNW DAVADQFSL S LRTLHRQLKQ QTGLTPQRYL NRLRLMKARH LL

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--