

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Escherichia coli NADH-quinone oxidoreductase subunit C/D(nuoC), Escherichia Coli Preis auf Anfrage CSB-EP482353ENM

Artikelname	Recombinant Escherichia coli NADH-quinone oxidoreductase subunit C/D(nuoC), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP482353ENM
Hersteller Artikelnummer	CSB-EP482353ENM
Alternativnummer	CSB-EP482353ENM-1, CSB-EP482353ENM-100, CSB-EP482353ENM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	E. coli
Tag	The tag type will be determined during production process.
UniProt	B7LAU8
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MVNNMTDLTA QEPAWQTRDH LDDPVIGELR NRFGPDAFTV QATRTGVPVV WIKREQLLEV GDFLKKLPKP YVMLFDLHGM DERLRTHREG LPAADFSVFY HLISIDRNRD IMLKVALAEN DLHVPTFTKL FPNANWYERE TWDLFGITFD GHPNLRIMM PQTWKGHPLR KDYPARATEF SPFELTKAKQ DLEMEALTFK PEEWGMKRG T ENEDFMFLNL GPNHPSAHGA FR

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--