

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Pseudomonas aeruginosa* Exodeoxyribonuclease 7 large subunit(xseA),*Escherichia Coli* Preis auf Anfrage CSB-EP483949EZY

Artikelname	Recombinant <i>Pseudomonas aeruginosa</i> Exodeoxyribonuclease 7 large subunit(xseA), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP483949EZY
Hersteller Artikelnummer	CSB-EP483949EZY
Alternativnummer	CSB-EP483949EZY-1, CSB-EP483949EZY-100, CSB-EP483949EZY-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B7UXJ1
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MRNDPFQRLG LDREVLTVSQ LNQRARLLLE DVFPQVWVEG ELSNLARPAS GHVYFTLKDS NAQIRCALFR QNALRVRQAL RDGLAVKVRG KISLFEGRGD YQLIADTVEP AGDGALRLAF EALKEKLAGE GLFASERKRL LPAHPRRIGI VSSPSGAVIR DIISVFRRRA PQVELTLVPT AVQGREAVAQ IVRALQLADR QGFDALILAR GGSLEDLWC FNEEAVARAV AA

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--