

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Streptomyces clavuligerus Insertion element IS116 uncharacterized 44.8 kDa protein, Escherichia Coli Preis auf Anfrage CSB-EP324855FOA

Artikelname	Recombinant Streptomyces clavuligerus Insertion element IS116 uncharacterized 44.8 kDa protein, Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP324855FOA
Hersteller Artikelnummer	CSB-EP324855FOA
Alternativnummer	CSB-EP324855FOA-1, CSB-EP324855FOA-100, CSB-EP324855FOA-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P19834
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSTRHDRIWV GIDAGKGHHW AVAVDADGET LFSTKVINDE AQLVTLIETA REREEVRWAV DISGRASTLL LALLVAHGQN VVYVPGRTVN RMSGAYKGEG KTDAKDARVI ADQARMRRDF APLDRPPELV TTLRLLTNHR ADLIADRVRL INLRDLLTG ICPALERAFF YSAAKGPVVM LNEYQTPAAL RRTGVKRLTT WLGRRKVRDA DTVAAKAIEA ARTQQVVLPG EK

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--