

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Vibrio splendidus* UDP-2,3-diacylglycosamine hydrolase(IpxH),*Escherichia Coli* Preis auf Anfrage CSB-EP492760VFH

Artikelname	Recombinant <i>Vibrio splendidus</i> UDP-2,3-diacylglycosamine hydrolase(IpxH), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP492760VFH
Hersteller Artikelnummer	CSB-EP492760VFH
Alternativnummer	CSB-EP492760VFH-1, CSB-EP492760VFH-100, CSB-EP492760VFH-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	B7VGZ3
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MKTYFISDLH LAPSRQDITD CFLTFMKNEA LEADALYVLG DLFEFWIGDD DTSEFANSIR QTFIDLVNTG VPCYFTQGNR DFLVGKKFAK QTGVQLLDEV STIDIYGQKA VVLHGDTLCT EDIKYLAFRE KVHQPWLQWV FNRIPFFIKK KIVSKVQSDI KDDKQTKSLD IMDVTQQEVE NVMKQHNVDL MIHGHTHRPN IHSFSANNCT KTRIVLGDWY TQGSVLVFTP QS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--