

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Paracoccus denitrificans Protein moxJ(moxJ),Escherichia Coli Preis auf Anfrage CSB-EP341105PBZ

Artikelname	Recombinant Paracoccus denitrificans Protein moxJ(moxJ),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP341105PBZ
Hersteller Artikelnummer	CSB-EP341105PBZ
Alternativnummer	CSB-EP341105PBZ-1, CSB-EP341105PBZ-100, CSB-EP341105PBZ-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	P29900
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	DTTNLR VCASTKDA PF SDAQGAGFEN KIAQVLADEM GATLDLVMLE KDAIYLVRDG IEKDLCDVLV GVDAGDERLL TTRPYRSGY AFVTRQDRNF EGDKWQDVDQ EGFDTFSYRL HSPAETILKY TGRYEYNLIY QASLTNFEDR RNKYTQVEAS RVITEVADGG ADLAIVFAPE AARYVRDSRE PLRMTLITNE IERSDGVII LQYSQSVGVS KTHPELLGPI EQALQS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--