

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Clostridium botulinum Holliday junction ATP-dependent DNA helicase RuvB(ruvB),Escherichia Coli Preis auf Anfrage CSB-EP501209DUE

Artikelname	Recombinant Clostridium botulinum Holliday junction ATP-dependent DNA helicase RuvB(ruvB),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP501209DUE
Hersteller Artikelnummer	CSB-EP501209DUE
Alternativnummer	CSB-EP501209DUE-1, CSB-EP501209DUE-100, CSB-EP501209DUE-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Clostridium
Tag	The tag type will be determined during production process.
UniProt	C3KTD2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MENRMVTPFD VEDDREQYSL RPTTLKEYIG QKKVKANLDI FIKAAKKRSE SLDHVLFYGP PGLGKTTLAN IANEMTGNL KVTSGPAIEK AGDLAAILTS LTDYDVLFDI EIHLNRSIE EILYPAMEDY ALDIVIGKGA AAKSIRLDLP KFTLIGATTR VGLLTSPLRD RFGMLCAMEF YTDEELMEIV VRSAAILNVN ICREAAFEIG KRSRGTPRIA NRLLKRVRDY CD

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--