

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Hamiltonella defensa* subsp. *Acyrtosiphon pisum* Outer-membrane lipoprotein lolB(lolB),*Escherichia Coli* Preis auf Anfrage CSB-EP504382HUG

Artikelname	Recombinant <i>Hamiltonella defensa</i> subsp. <i>Acyrtosiphon pisum</i> Outer-membrane lipoprotein lolB(lolB), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP504382HUG
Hersteller Artikelnummer	CSB-EP504382HUG
Alternativnummer	CSB-EP504382HUG-1, CSB-EP504382HUG-100, CSB-EP504382HUG-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	C4K7Y9
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	CSALKAPES SSVLKNHIAS DEWQEYQHQL KQIQQFQIQG SVAYFSDEKK AYARFFWQQY SPKNYHLLLL SPLGQTEFEL KVTNGRVDMA KYKDQGEIKG DAEEILFKLT GIPILEHLS RWIVGASSDA DEILNRQSR LKTLIHHKKE QTWKVYYQAY NTKITPILPE RLELYLISPN HQDQRMTLKI NHWILK

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--