

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Desulfovibrio magneticus* DNA repair protein recO(recO), *Escherichia Coli* Preis auf Anfrage CSB-EP505444DJP

Artikelname	Recombinant <i>Desulfovibrio magneticus</i> DNA repair protein recO(recO), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP505444DJP
Hersteller Artikelnummer	CSB-EP505444DJP
Alternativnummer	CSB-EP505444DJP-1, CSB-EP505444DJP-100, CSB-EP505444DJP-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	C4XS52
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MEYSDQALIL AVKRFREIDA WVRLLSPTRG VYTAFAFGGL KSRRRFLGCL DPLNHVQFKV RRSGYRGYHC LAEGKLLDAP RQLRSHQPRL GMAVNCLKFF EAAPVAPGSF AQAYGLMRAM LSTLDAADep SPLFPLLFR RMTFLHGMLP ACGQCAVCGQ PLGHDGAVCH VEEGRVACPD CRAAASGGVH ARLGGEALAL LASAVEQGPE QWAACRPHPA AGREFSRAVD LL

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--