

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Thermococcus sibiricus Ribonuclease Z(rnz), Escherichia Coli Preis auf Anfrage CSB-EP513666TMT

Artikelname	Recombinant Thermococcus sibiricus Ribonuclease Z(rnz), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP513666TMT
Hersteller Artikelnummer	CSB-EP513666TMT
Alternativnummer	CSB-EP513666TMT-1, CSB-EP513666TMT-100, CSB-EP513666TMT-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	C6A2P3
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MLEVIFLGTG GIMPTKERNV PAIALKYGGE IILWDIGEGT LRQLDTTKLS PMKIEKIFIT HFHGDHYLGL MSLIQTMTLW KRERPLYIYG PKYTFEFIQN YLKSGFFRPS FEIHIHELGE TRLKFKDYEI WSKVDHGIP ALGYVFKEND KRGNFDVEKI RELGLKPGPW MKKLEIKGKI EINGKLIYLE DVTGEPKKGV KIVYTGDETEQ CERIKLFSEADLLIHEATY LS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--