

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Shewanella sp. DNA polymerase IV(dinB), Escherichia Coli Preis auf Anfrage CSB-EP367286STS

Artikelname	Recombinant Shewanella sp. DNA polymerase IV(dinB), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP367286STS
Hersteller Artikelnummer	CSB-EP367286STS
Alternativnummer	CSB-EP367286STS-1, CSB-EP367286STS-100, CSB-EP367286STS-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A0KTR5
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MRKIIHIDMD CYFAAVEMRD FPEYRGKPLA VGGSRVQRGV ISTCNYEARK FGVRSAMATG YALKLCPDLI LVPGRMQVYK EVSQQIRAIF CRYTELIEPL SLDEAYLDVS DCKLFGKSAT LIAEAIRRD I LAETGLTASA GVAPIKFLAK VASDLNKPNG QCVIPPDEVA EFVKSLSLRK IPGVGKVTA E KLSSLGLNTC ADVQAYPKQE LIARFGKFGT VLVERAHGID ER

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--