

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Aquifex aeolicus Glycerol-3-phosphate dehydrogenase [NAD(P)+], Escherichia Coli Preis auf Anfrage CSB-EP525468DNV

Artikelname	Recombinant Aquifex aeolicus Glycerol-3-phosphate dehydrogenase [NAD(P)+], Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP525468DNV
Hersteller Artikelnummer	CSB-EP525468DNV
Alternativnummer	CSB-EP525468DNV-1, CSB-EP525468DNV-100, CSB-EP525468DNV-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	O67555
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MHEIFYTGCG KVGNTLATHL DRKGHEVLLY DINEEVVKKI NEGHPYVEG IRFKNVKATT DLNQINEFKN LICALPVQVI PKVITKVNLIK GKNFISASKG IIHENLKRIS QLVKEIEPKL KFFVLSGPSF AEEVSKGLPT AIVLAYEDKE EAMKLQDALD SENFNVLND DITGVELGGA LKNVIAIavg LSDGMGYGYN ARSAITRGL HEMANVGEKF GARRETFEGL SG

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--