

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant *Aspergillus oryzae* Endo-1,4-beta-xylanase F1(xynF1),*Escherichia Coli* Preis auf Anfrage CSB-EP526074DPA

|                          |                                                                                                                                                                                                                                                                                |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant <i>Aspergillus oryzae</i> Endo-1,4-beta-xylanase<br>F1(xynF1), <i>Escherichia Coli</i> Preis auf Anfrage                                                                                                                                                           |
| Artikelnummer            | CSB-EP526074DPA                                                                                                                                                                                                                                                                |
| Hersteller Artikelnummer | CSB-EP526074DPA                                                                                                                                                                                                                                                                |
| Alternativnummer         | CSB-EP526074DPA-1, CSB-EP526074DPA-100, CSB-EP526074DPA-20                                                                                                                                                                                                                     |
| Hersteller               | Cusabio                                                                                                                                                                                                                                                                        |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                               |
| Spezies Reaktivität      | Bacteria                                                                                                                                                                                                                                                                       |
| Tag                      | The tag type will be determined during production process.                                                                                                                                                                                                                     |
| UniProt                  | <a href="#">O94163</a>                                                                                                                                                                                                                                                         |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                    |
| Quelle                   | E.coli                                                                                                                                                                                                                                                                         |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                    |
| Formulierung             | Lyophilized powder                                                                                                                                                                                                                                                             |
| Sequenz                  | A VISRQAAASI NDAFVAHGKK YFGTCSDQAL LQNSQNEAIV<br>RADFGQLTPE NSMKWDALEP SQGSFSFAGA DFLADYAKTN<br>NKLVRGHTLV WHSQLPSWVQ GITDKDTLTE VIKNHITTIM QRYKGQIYAW<br>DVVNEIFDED GTLRDSVFSQ VLGEDFVRIA FETAREADPN AKLYINDYNL<br>DSADYAKTKG MVS YVKKWLD AGVPIDGIGS QSHYSANGFP<br>VSGAKGALTA |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|