

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Phaeodactylum tricornutum Photosystem I iron-sulfur center(psaC),Escherichia Coli Preis auf Anfrage CSB-EP371556EUF

Artikelname	Recombinant Phaeodactylum tricornutum Photosystem I iron-sulfur center(psaC),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP371556EUF
Hersteller Artikelnummer	CSB-EP371556EUF
Alternativnummer	CSB-EP371556EUF-1, CSB-EP371556EUF-100, CSB-EP371556EUF-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A0T0L2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	SHTVKIYDT CIGCTQCVRA CPTDVLEMVP WDGCKAGQIA SSPRVEDCVG CKRCETACPT DFLSVRVYLG AETTRSLGLA Y

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.