

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Escherichia coli O1:K1 / APEC HTH-type transcriptional activator AIIIS(aIIIS), Escherichia Coli Preis auf Anfrage CSB-EP371728EJE

Artikelname	Recombinant Escherichia coli O1:K1 / APEC HTH-type transcriptional activator AIIIS(aIIIS), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP371728EJE
Hersteller Artikelnummer	CSB-EP371728EJE
Alternativnummer	CSB-EP371728EJE-1, CSB-EP371728EJE-100, CSB-EP371728EJE-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	E. coli
Tag	The tag type will be determined during production process.
UniProt	A1A8H0
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MFDPETLRTF IAVAETGSFS KAAERLCKTT ATISYRIKLL EENTGVALFF RTTRSVTLTA AGEHLLCQAR DWLGWLESMP SELQQVNDGV ERQVNIVINN LLYNPQAVAR LLAWLNERYP FTQFHISRQI YMGVWDSLLY EGFSLAIGVT GTEALANTFS LDPLGSVQWR FVMAADHPLA NVEEPLTEAQ LRRFPAVNIE DSARTLTKRV AWRLPGQKEI IVPDMETKIA AH

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--