

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Caulobacter sp. ATP-dependent protease ATPase subunit HslU(hslU),Escherichia Coli Preis auf Anfrage CSB-EP536190DQR

Artikelname	Recombinant Caulobacter sp. ATP-dependent protease ATPase subunit HslU(hslU),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP536190DQR
Hersteller Artikelnummer	CSB-EP536190DQR
Alternativnummer	CSB-EP536190DQR-1, CSB-EP536190DQR-100, CSB-EP536190DQR-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	B0T7A8
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTEFSPREIV SELDRYIVGH TEAKKAVAVA LRNRWRRRRV PADLRDEVTP KNILLIGPTG VGKTEIARRL AKLAQAPFLK VEATKFTEVG YVGRDVDQIV RDLVESALSM VRDKRRGGVK AKAEAAAAEER ILDALTGPGS TAARESFRKK LRAGELDDKE VELQLADTGG GGMFEIPGQP GASVLNLSEM MKSFGGGRTK THKTTVVGAW APLIAEESDK LLDQEALTQE AL

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--