

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Francisella tularensis subsp. tularensis ATP-dependent Clp protease proteolytic subunit(clpP),Escherichia Coli Preis auf Anfrage CSB-EP392120FDX

Artikelname	Recombinant Francisella tularensis subsp. tularensis ATP-dependent Clp protease proteolytic subunit(clpP),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP392120FDX
Hersteller Artikelnummer	CSB-EP392120FDX
Alternativnummer	CSB-EP392120FDX-1, CSB-EP392120FDX-100, CSB-EP392120FDX-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A4IYB6
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MITNNLVPTV IEKTAGGERA FDIYSRLLKE RIVFLNGEVN DHSANLVIAQ LLFLESEDPD KDIYFYINSP GGMVTAGMGV YDTMQFIKPD VSTICIGLAA SMGSLLLAGG AKGKRYSLPS SQIMIHQPLG GFRGQASDIE IHAKNILRIK DRLNKVLAHH TGQDLETIVK DTDNRDNFMMMA DEAKAYGLID HVIESREAI K

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--