

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Mycobacterium abscessus Non-canonical purine NTP pyrophosphatase (MAB_1485), Escherichia Coli Preis auf Anfrage CSB-EP544522MVF

Artikelname	Recombinant Mycobacterium abscessus Non-canonical purine NTP pyrophosphatase (MAB_1485), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP544522MVF
Hersteller Artikelnummer	CSB-EP544522MVF
Alternativnummer	CSB-EP544522MVF-1, CSB-EP544522MVF-100, CSB-EP544522MVF-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	B1MLZ4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MKILVASRNP KKLAELSRVL ESSGVSGVEL VSLTDVPEYE EVPETGASFE DNALIKAREG VKHTGLACVA DDSGLAVDAL NWMPGVLSAR WSGRHGDDAA NTALLLAQLS DIPDERRGAA FVSACALVTP EGEEVVVEGR WKGSIARIPA GQNGFGYDPI FVPRGGLRTA AELTPEEKDA VSHRGRALAA LLPMLRNLVN LGRTAP

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.