

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Recombinant Clavibacter michiganensis subsp. michiganensis Imidazoleglycerol-phosphate dehydratase(hisB),Escherichia Coli Preis auf Anfrage CSB-EP397112CWT**

|                          |                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Clavibacter michiganensis subsp. michiganensis<br>Imidazoleglycerol-phosphate dehydratase(hisB),Escherichia Coli Preis<br>auf Anfrage |
| Artikelnummer            | CSB-EP397112CWT                                                                                                                                   |
| Hersteller Artikelnummer | CSB-EP397112CWT                                                                                                                                   |
| Alternativnummer         | CSB-EP397112CWT-1, CSB-EP397112CWT-100, CSB-EP397112CWT-<br>20                                                                                    |
| Hersteller               | Cusabio                                                                                                                                           |
| Kategorie                | Proteine/Peptide                                                                                                                                  |
| Spezies Reaktivität      | Bacteria                                                                                                                                          |
| Tag                      | The tag type will be determined during production process.                                                                                        |
| UniProt                  | <a href="#">A5CSK8</a>                                                                                                                            |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                       |
| Quelle                   | E.coli                                                                                                                                            |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                       |
| Formulierung             | Lyophilized powder                                                                                                                                |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                             |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz                | MSTLARTAHV TRQTSESTID LQLDLDTGA SEISTSVPFY DHMLTAFKHS<br>SLTDLRVTAT GDTHIDVHHT VEDVGIVLGQ AIREALGDKS GIARFGDALV<br>PLDEALVQSV VDISGRPFLV HSGEPAGFEM HLIGGHFTGS MVRHVFEAIT<br>FHAGLTVHVT VLGGRDPHHI AEAEFKSFAR AFRQAKELDP RVSGIPSTKG<br>AL                                                                                                                                                                   |
| Anwendungsbeschreibung | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |