

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Vibrio cholerae serotype O1 UDP-N-acetylmuramoylalanine--D-glutamate ligase(murD),Escherichia Coli Preis auf Anfrage CSB-EP399670VEY

Artikelname	Recombinant Vibrio cholerae serotype O1 UDP-N-acetylmuramoylalanine--D-glutamate ligase(murD),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP399670VEY
Hersteller Artikelnummer	CSB-EP399670VEY
Alternativnummer	CSB-EP399670VEY-1, CSB-EP399670VEY-100, CSB-EP399670VEY-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5F5N1
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MDRWQGIQHV VVGLGITGL SVVNYLRKYH PSVTVQVIDT REIPPGQEQL SSDVALHRGG WNLEWLLNAD LVVTNPGIAL ATPEIQQVLA AGIPVVGDI LFAWHVDTPV IAITGSNGKS TVTDLSGVLA NAAGVKA AVG GNIGVPALDL ISPDVELYVL ELSSFQLETT SSLKLKAAAF LNLSEDHMDR YQGMGDYRQA KLRIFDHAET AVVNADDTQT FPDHAAHLQV VT

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.