

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Flavobacterium johnsoniae NADH-quinone oxidoreductase subunit D(nuoD), Escherichia Coli Preis auf Anfrage CSB-EP400017FDT

Artikelname	Recombinant Flavobacterium johnsoniae NADH-quinone oxidoreductase subunit D(nuoD), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP400017FDT
Hersteller Artikelnummer	CSB-EP400017FDT
Alternativnummer	CSB-EP400017FDT-1, CSB-EP400017FDT-100, CSB-EP400017FDT-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	A5FKI7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSELLLPPEH RYAKIIKERH NEDGSELSVL NLGPTHPTH GIFQNILLMD GERILEAEPT IGYIHRAFEK IAENRPFYQI TPLTDRMNYC SSPINMGWW MTLEKLLGIE VPKRAQYLRV IVMEALARITD HIICNSILGV DTGAYTGFLY VFQFREKIYE IYEEICGARL TTNMGRIGGF ERDWSPEAFR KLDKFLEEF PAWKEFENLF ERNRIFLDRT VNVGGITAEK AM

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.