

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Alcanivorax borkumensis Methionyl-tRNA formyltransferase(fmt),Escherichia Coli Preis auf Anfrage CSB-EP611807AAAM

Artikelname	Recombinant Alcanivorax borkumensis Methionyl-tRNA formyltransferase(fmt),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP611807AAAM
Hersteller Artikelnummer	CSB-EP611807AAAM
Alternativnummer	CSB-EP611807AAAM-1, CSB-EP611807AAAM-100, CSB- EP611807AAAM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q0VTE2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MKPLRLAFAG TPDFAAASLQ AVLDNGHQVV AVLTQPDRAA GRGKKLQQSP VKQLAHSQGI TVLQPENLKG EAIHQQLRDL NLDALVVVAY GLIIPQAVLD MPRLGCLNVH GSLLPRWRGA APIQRAITAG DTETGNTIMQ MEAGLDTGPM LLSLESLPIGD SETGGELHDR LAAQGARLLV TVLQDLEKYL ANATQPDEG ITYAHKLSKA EARLDFRLPT GALYNRIRAF NP

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.