

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Haemophilus influenzae Membrane-bound lytic murein transglycosylase F(mltF), Escherichia Coli Preis auf Anfrage CSB-EP404075HSE

Artikelname	Recombinant Haemophilus influenzae Membrane-bound lytic murein transglycosylase F(mltF), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP404075HSE
Hersteller Artikelnummer	CSB-EP404075HSE
Alternativnummer	CSB-EP404075HSE-1, CSB-EP404075HSE-100, CSB-EP404075HSE-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A5UG36
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	ID MVFPWQFLRH TEENHYTAIQ ARGSLYVGTI NNQISYFINK DSERGFYEYEL AKAFADTLGV ELEMKIFDNQ EQLFDELNKH NIDLAAAHLL YHPKNAERFQ IGPAYHSASW QLAYRKNENR PKNLGNVKKD IYISNNLALE ETLKEKQKQY PQLIWKRNQA LTQEELLQL AEGKIPYVIA NSIDIAAIQQ IKPELAIAFD ITDEANVHWY LPNSPYRDLQ TALLNFMNNA

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--