

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Haemophilus influenzae UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase(murE),Escherichia Coli Preis auf Anfrage **CSB-EP404923HSE**

Artikelname	Recombinant Haemophilus influenzae UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase(murE),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP404923HSE
Hersteller Artikelnummer	CSB-EP404923HSE
Alternativnummer	CSB-EP404923HSE-1, CSB-EP404923HSE-100, CSB-EP404923HSE-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A5UIQ7
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MKKLTALFNL PELKNDIELH NMVLDSRKVK AGDLFVAIKG HQVDGNQFID SALHSGASAV VSETELSSEH LTVEFIGNVP VVKYHQLARH LSSLADVFD SPSKNLTLVG VTGTNGKTTI SQLLAQWAEI LGHRAAVMGT IGNGLFRQIV EAKNTTGSAV EIQSSLSTFK HAGADFTSIE VSSHGLAQHR VEALHFKAAL FTNLTRDHLD YHQSMENYAA AKKRLFTELD TQ
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.