

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant Yersinia pestis bv. Antiqua HTH-type transcriptional activator rhaR(rhaR),Escherichia Coli Preis auf Anfrage CSB-EP618471YAF

|                          |                                                                                                                                                                                                                                                                                    |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Yersinia pestis bv. Antiqua HTH-type transcriptional activator rhaR(rhaR),Escherichia Coli Preis auf Anfrage                                                                                                                                                           |
| Artikelnummer            | CSB-EP618471YAF                                                                                                                                                                                                                                                                    |
| Hersteller Artikelnummer | CSB-EP618471YAF                                                                                                                                                                                                                                                                    |
| Alternativnummer         | CSB-EP618471YAF-1, CSB-EP618471YAF-100, CSB-EP618471YAF-20                                                                                                                                                                                                                         |
| Hersteller               | Cusabio                                                                                                                                                                                                                                                                            |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Bacteria                                                                                                                                                                                                                                                                           |
| Tag                      | The tag type will be determined during production process.                                                                                                                                                                                                                         |
| UniProt                  | <a href="#">Q1C0W2</a>                                                                                                                                                                                                                                                             |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                        |
| Quelle                   | E.coli                                                                                                                                                                                                                                                                             |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                        |
| Formulierung             | Lyophilized powder                                                                                                                                                                                                                                                                 |
| Sequenz                  | MRAPLLLES R DYLLSEQMPV AVTNRY PQET FVEH THQFCE IVIVWRGNGL<br>HVLNDHPYRI TCGDV FYIQA ADHHSY ES VH DLVLDNIIYC PERLHLNAQW<br>HKLLPPLGPE QNQGYWRLTT QGMAQARPII QQLAQESRKT DSWSIQLTEV<br>LLLQLAIVLK RHR YRAEQAH LLPDGEQLDL IMSALQQSLG AYFDMADFCH<br>KNQLVERSLK QLFRQQTGMS ISHYLRQIRL CH |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|