

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Shewanella baltica Peptidase E(pepE), Escherichia Coli Preis auf Anfrage CSB-EP408051STM

Artikelname	Recombinant Shewanella baltica Peptidase E(pepE), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP408051STM
Hersteller Artikelnummer	CSB-EP408051STM
Alternativnummer	CSB-EP408051STM-1, CSB-EP408051STM-100, CSB-EP408051STM-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6WPJ2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTINALLSS SRVGDTPLYA HAIPFIKPLT TNAQKWIFIP YAGVSMYSYN YLASVVTGLS ELELDISGIH QHPDPQQAIAK DADGILIGGG NTFHLLHQLY RYDLVTLIGE QVALGKPYIG WSAGSNVSGQ SIRTNDMPPI IEPPSFKALN LLPFQLNPHY SNYQAPGHNG ETRAQRLLEF TKVDPLTPVV GIVEGSALWR QGDKLSLLGD QPAYLFCGEQ QEIPVVGSD LS

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--