

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Burkholderia cenocepacia ATP-dependent dethiobiotin synthetase BioD(bioD), Escherichia Coli Preis auf Anfrage CSB-EP622302BAAG

Artikelname	Recombinant Burkholderia cenocepacia ATP-dependent dethiobiotin synthetase BioD(bioD), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP622302BAAG
Hersteller Artikelnummer	CSB-EP622302BAAG
Alternativnummer	CSB-EP622302BAAG-1, CSB-EP622302BAAG-100, CSB-EP622302BAAG-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q1BT35
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTAPLSVFVT GTDTEIGKTF VSAAMLHGFA RHGLRAAALK PVAAGAYERD GVWRNEDADQ LDAAANVALP PELRTPFLK APAAPHIVAA QEGVTLDLDT IVACHREALT RADVVVVEGV GGFRVPLNDT QDTADLAVAL GLPVVLVVG RLGCISHALL TADAIRQRGL TLAGWVANHV DPAMSFADEN VATIRDWLAR EHRAPLIGRI AHMTPAAPES AAAMLIDIAAL VE

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--