

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Klebsiella pneumoniae subsp. pneumoniae (3R)-hydroxymyristoyl-[acyl-carrier-protein] dehydratase, Escherichia Coli Preis auf Anfrage CSB-EP408445KAX

Artikelname	Recombinant Klebsiella pneumoniae subsp. pneumoniae (3R)-hydroxymyristoyl-[acyl-carrier-protein] dehydratase, Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP408445KAX
Hersteller Artikelnummer	CSB-EP408445KAX
Alternativnummer	CSB-EP408445KAX-1, CSB-EP408445KAX-100, CSB-EP408445KAX-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6T4Y2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTTDTHTLHI EEILELLPHR YPFLLVDRV LDFEEGRFLRA VKNVSVNEPF FQGHFPGKPI LPGVLILEAM AQATGILAFK SVGKLEPGEL YYFAGIDEAR FKRPVVPGDQ MIMVTFEKT RRLTRFKGV ALVDGKVVCE ATMMCAR SRE A

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--