

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Klebsiella pneumoniae subsp. pneumoniae Ribosomal protein S12 methylthiotransferase RimO(rimO), Escherichia Coli Preis auf Anfrage CSB-EP408453KAX

Artikelname	Recombinant Klebsiella pneumoniae subsp. pneumoniae Ribosomal protein S12 methylthiotransferase RimO(rimO), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP408453KAX
Hersteller Artikelnummer	CSB-EP408453KAX
Alternativnummer	CSB-EP408453KAX-1, CSB-EP408453KAX-100, CSB-EP408453KAX-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6T6T1
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MSNVTHQPKI GFVSLGCPKN LVDSERILTE LRTEGYDVVP TYDNADMVIV NTCGFIDSAV QESLEAIGEA LKENGKVIVT GCLGAKEDQI REVHPKVLEI TGPHSYEQVL EHVHHYTPKP KHNPFLSLVP EQGVKLTPRH YAYLKISEGC NHRCTFCIIP SMRGDLVSRP IGEVLAEAKR LADAGVKELL VISQDTSAYG VDVKHRTGFH NGMPVKTSMV SLCEELAKLG IW
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.