

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Pseudomonas aeruginosa* Succinyl-diaminopimelate desuccinylase(dapE),*Escherichia Coli* Preis auf Anfrage CSB-EP409452PZL

Artikelname	Recombinant <i>Pseudomonas aeruginosa</i> Succinyl-diaminopimelate desuccinylase(dapE), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP409452PZL
Hersteller Artikelnummer	CSB-EP409452PZL
Alternativnummer	CSB-EP409452PZL-1, CSB-EP409452PZL-100, CSB-EP409452PZL-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	A6V936
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MTASSPSLSPTLELACELIRPSVTPVDADCQALMMRRLEAAGFALEPMRIEEVDNFWARRGGDGPVLCFAGHTDVVPTGPVQAWQHQPFDALIDEQGMLCGRGAADMKGSLASMIVAVERFVADHPGHKGAI AFLITSD EEGPAHHGTKAVVERLAARGERLDWCIVGEPSS TSLVGDVVKNGRRGSLGARLTIRGVQGHVAYPHLAKNPIHLAAPALAEAAEHWDDGNA

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--