

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Recombinant *Klebsiella pneumoniae* subsp. *pneumoniae* Non-specific ribonucleoside hydrolase rihC(rihC), *Escherichia Coli* Preis auf Anfrage CSB-EP411642KAX

|                          |                                                                                                                                                                                                                                                                             |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant <i>Klebsiella pneumoniae</i> subsp. <i>pneumoniae</i> Non-specific ribonucleoside hydrolase rihC(rihC), <i>Escherichia Coli</i> Preis auf Anfrage                                                                                                               |
| Artikelnummer            | CSB-EP411642KAX                                                                                                                                                                                                                                                             |
| Hersteller Artikelnummer | CSB-EP411642KAX                                                                                                                                                                                                                                                             |
| Alternativnummer         | CSB-EP411642KAX-1, CSB-EP411642KAX-100, CSB-EP411642KAX-20                                                                                                                                                                                                                  |
| Hersteller               | Cusabio                                                                                                                                                                                                                                                                     |
| Kategorie                | Proteine/Peptide                                                                                                                                                                                                                                                            |
| Spezies Reaktivität      | Bacteria                                                                                                                                                                                                                                                                    |
| Tag                      | The tag type will be determined during production process.                                                                                                                                                                                                                  |
| UniProt                  | <a href="#">A6T4G4</a>                                                                                                                                                                                                                                                      |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                                                                                                                                                                                 |
| Quelle                   | E.coli                                                                                                                                                                                                                                                                      |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                 |
| Formulierung             | Lyophilized powder                                                                                                                                                                                                                                                          |
| Sequenz                  | MRLPIILDTD PGIDDAAAIA AALFAPELDL QLMTTVAGNV SVEKTTRNAL<br>QLLHFWNADV PLAQGASMPV VRPLRDAASV HGESGMEGYD<br>FVEHQEQPLA KPAFQAIRDA LMHAAEPITL VAIGPLTNIA LLLTQYPECV<br>FNIRRLVIMG GSAGRGNFTP NAEFNIAIDP EAAAKVFHSG LEIVMCGLDV<br>TNRALLAADY LATLPTLNQT GKMLHALFSH YRSGSMSSGL RM |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|