

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Xanthomonas campestris pv. campestris UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase(murE),Escherichia Coli Preis auf Anfrage CSB-EP819008XAY

Artikelname	Recombinant Xanthomonas campestris pv. campestris UDP-N-acetylmuramoyl-L-alanyl-D-glutamate--2,6-diaminopimelate ligase(murE),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP819008XAY
Hersteller Artikelnummer	CSB-EP819008XAY
Alternativnummer	CSB-EP819008XAY-1, CSB-EP819008XAY-100, CSB-EP819008XAY-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q8PCK4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSRAMALSQL LPDVALARDV QVSGLVMDSR AVRPGDAFVA IAGFGAHGLG FAAQALASGA SAILFEPPAP ADLPVPADAI AVPGLRARLG VLADHFH GAP SQAMRMIGVT GTNGKTSTVQ LLAQALTLLG TPTGTIGTLG VGLYGAAVPT GFTTPLVLQT HAELAQLRDA GAQAVAMEVS SHALDQGRVD AVQFDVAVFT NLTRDHLDYH GDMAQYGA AK AKLFARAGLR AA

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--