

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Ehrlichia chaffeensis Lipoprotein-releasing system ATP-binding protein LoID(loID),Escherichia Coli Preis auf Anfrage CSB-EP646195EJS

Artikelname	Recombinant Ehrlichia chaffeensis Lipoprotein-releasing system ATP-binding protein LoID(loID),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP646195EJS
Hersteller Artikelnummer	CSB-EP646195EJS
Alternativnummer	CSB-EP646195EJS-1, CSB-EP646195EJS-100, CSB-EP646195EJS-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q2GHT4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSTVFALANI YKSFGRNKEV PVIVNATLQI KKGEIVALIG KSGSGKSTLL HIAGLLDTPS SGSILLNIE CTDKTSKDK TYLRRHFLGF IYQFHHLLQE FSVLENVMLP QIITGKSKI AKKNAMEILE RVRLQDKLSM PISQLSGGER QRVAIARSLI NCPLIVLADE PTGSLDNNTA LEVFSLLQEY AKEKNIAILL ATHNYSLAQK ACRIVKIDSG ILRSYSIDES GF

Anwendungsbeschreibung

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.