

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Shigella boydii serotype 4 Probable L-ascorbate-6-phosphate lactonase ulaG(ulaG),Escherichia Coli Preis auf Anfrage CSB-EP656052SGE

Artikelname	Recombinant Shigella boydii serotype 4 Probable L-ascorbate-6-phosphate lactonase ulaG(ulaG),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP656052SGE
Hersteller Artikelnummer	CSB-EP656052SGE
Alternativnummer	CSB-EP656052SGE-1, CSB-EP656052SGE-100, CSB-EP656052SGE-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q31TC3
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MSKVKSITRE SWILSTFPEW GSWLNEEIEQ EQVAPGTFAM WWLGCTGIWL KSEGGTNVCV DFWCGTGKQS HGNPLMKQGH QMORMAGVKK LQPNLRTPPF VLDPF AIRQI DAVLATHDHN DHIDVNVA VMQNCADDVP FIGPKTCVDL WIGWGVPKER CIVVKPGDVV KVKDIEIHAL DAFDRTALIT LPADQKAAGV LPDGMDDRAV NYLFKTPGGS LYHSGDSHYS NY

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--