

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Methanopyrus kandleri 5-formaminoimidazole-4-carboxamide-1-(beta)-D-ribofuranosyl 5-monophosphate synthetase(purP),Escherichia Coli Preis auf Anfrage CSB-EP855178MSI

Artikelname	Recombinant Methanopyrus kandleri 5-formaminoimidazole-4-carboxamide-1-(beta)-D-ribofuranosyl 5-monophosphate synthetase(purP),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP855178MSI
Hersteller Artikelnummer	CSB-EP855178MSI
Alternativnummer	CSB-EP855178MSI-1, CSB-EP855178MSI-100, CSB-EP855178MSI-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q8TXS2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MLSREEILEI LSEYDLEDLS IATLGSH TAL HILKGAKEEG LRSVVVCEEG RTTPYERLGV ADEIIVVERF QDMLDEEVQE RLRELNAIVV PHGSFVAYVG LDGIENEFCV PMFGNRRLLR WESERSLERK LLKRA GVKVP KVF DSPEDID RPVIVKFPGA RGGRGYFICS DPEEFEEKAE RLIEDGVIDE EDLEQA HIEE YVVG TNFCVH YFRSVVEDTV EVLGMDRRYE TN
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.