

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Streptococcus pyogenes serotype M28 Glycerol-3-phosphate dehydrogenase [NAD(P)+](gpsA),Escherichia Coli Preis auf Anfrage CSB-EP670912SBAF

Artikelname	Recombinant Streptococcus pyogenes serotype M28 Glycerol-3-phosphate dehydrogenase [NAD(P)+](gpsA),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP670912SBAF
Hersteller Artikelnummer	CSB-EP670912SBAF
Alternativnummer	CSB-EP670912SBAF-1, CSB-EP670912SBAF-100, CSB-EP670912SBAF-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q48VF4
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MTKQKVAILG PGSWGTTALSQ VLNDNGHDVR LWGNIPDQIE EINTKHTNRH YFKDIVLDKN ITATLDLGQA LSDVDAVLV VPTKVTRLVA RQVAAILDHK VVMHASKGL EPETHERLST ILEEEIPAHF RSEVVVSGP SHAEETIVRD ITLITAASKD IEAAKYVQSL FSNHYFRLYT NTDVIGVETA GALKNIIAVG AGALHGLGYG DNAKAAVITR GLAEITRLGV KL
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.