

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Arabidopsis thaliana Peroxidase 56(PER56),Escherichia Coli Preis auf Anfrage CSB-EP873392DOA

Artikelname	Recombinant Arabidopsis thaliana Peroxidase 56(PER56),Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP873392DOA
Hersteller Artikelnummer	CSB-EP873392DOA
Alternativnummer	CSB-EP873392DOA-1, CSB-EP873392DOA-100, CSB-EP873392DOA-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Plant
Tag	The tag type will be determined during production process.
UniProt	Q9LXG3
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	LKVGFYDKA CPKAELIVKK SVFEAVKNDR TIAAPLLRMF FHDCFVRGCE GSVLLLELKNK KDEKNSIPNL TLRGFEIIDN VKAALEKECP GIVSCSDVLA LVARDAMVAL NGPSWEVETG RRDGLVTNIT EALLNLPSPF NNISSLITQF QSKGLDKKDL VVLSGGHTIG NGHCPQITNR LYNFTGKGDS DPNLDTEYAV KLRGKCKPTD TTTALEMDPG SFKTFDESYF KLV

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--