

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### **Recombinant Neisseria gonorrhoeae UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage CSB-EP682795NGA**

|                          |                                                                                                                          |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | Recombinant Neisseria gonorrhoeae UDP-N-acetylenolpyruvoylglucosamine reductase(murB),Escherichia Coli Preis auf Anfrage |
| Artikelnummer            | CSB-EP682795NGA                                                                                                          |
| Hersteller Artikelnummer | CSB-EP682795NGA                                                                                                          |
| Alternativnummer         | CSB-EP682795NGA-1, CSB-EP682795NGA-100, CSB-EP682795NGA-20                                                               |
| Hersteller               | Cusabio                                                                                                                  |
| Kategorie                | Proteine/Peptide                                                                                                         |
| Spezies Reaktivität      | Bacteria                                                                                                                 |
| Tag                      | The tag type will be determined during production process.                                                               |
| UniProt                  | <a href="#">Q5F9J9</a>                                                                                                   |
| Puffer                   | Tris/PBS-based buffer, 6% Trehalose, pH 8.0                                                                              |
| Quelle                   | E.coli                                                                                                                   |
| Reinheit                 | Greater than 85% as determined by SDS-PAGE.                                                                              |
| Formulierung             | Lyophilized powder                                                                                                       |

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                    |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sequenz                | MQPIRYRTDL TPYNTEFGLHA QARAFIALKH ADELRDIVRL PEFDRDRTLW<br>LGGGSNILLM QDYDGLVVHM ENKGIREIAR SDGMVLIEAQ AGEIWHDFVL<br>HTVALGLSGL ENSLIPGTV GASPVQNIGA YGVEAKDVIH SVRCFDLDE<br>TFVTLSNADC RFAYRESLFK QEGKGRYVIV SVVFALKTHF VPNLGYGDLA<br>AKVAELSAGR EATAKDVSDA VSAIRNSKLP DP                                                                                                                                         |
| Anwendungsbeschreibung | <p>We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.</p> |