

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant *Vibrio cholerae* serotype O1 UDP-3-O-[3-hydroxymyristoyl] N-acetylglucosamine deacetylase(IpxC),*Escherichia Coli* Preis auf Anfrage CSB-EP885190VEX

Artikelname	Recombinant <i>Vibrio cholerae</i> serotype O1 UDP-3-O-[3-hydroxymyristoyl] N-acetylglucosamine deacetylase(IpxC), <i>Escherichia Coli</i> Preis auf Anfrage
Artikelnummer	CSB-EP885190VEX
Hersteller Artikelnummer	CSB-EP885190VEX
Alternativnummer	CSB-EP885190VEX-1, CSB-EP885190VEX-100, CSB-EP885190VEX-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Spezies Reaktivität	Bacteria
Tag	The tag type will be determined during production process.
UniProt	Q9KPH2
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder
Sequenz	MIRQRTLKEI VKTTGVGLHS GRKVTLTLP AAANTGIIYR RTDVNPPVDF PADPASVRDT MLCTALVNDQ GVRISTVEHL NAALAGMGID NAIIEVDAPE IPIMDGSASP FVYLLQQAGI QTLNAPKRFI RIKKPVRIED GDKWAEFVPF NGFRMDFEIE FNHPAIDGDD QRLVDFDSSQ GFVKEISRAR TFGFMRDIEY LQSQNLCLGG SFDCAIVLDD YRILNEEGLR FD

Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
------------------------	--