

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant Bradyrhizobium japonicum Methylenetetrahydrofolate--tRNA-(uracil-5-)-methyltransferase TrmFO(trmFO), Escherichia Coli Preis auf Anfrage CSB-EP807215BVW

Artikelname	Recombinant Bradyrhizobium japonicum Methylenetetrahydrofolate--tRNA-(uracil-5-)-methyltransferase TrmFO(trmFO), Escherichia Coli Preis auf Anfrage
Artikelnummer	CSB-EP807215BVW
Hersteller Artikelnummer	CSB-EP807215BVW
Alternativnummer	CSB-EP807215BVW-1, CSB-EP807215BVW-100, CSB-EP807215BVW-20
Hersteller	Cusabio
Kategorie	Proteine/Peptide
Tag	The tag type will be determined during production process.
UniProt	Q89L21
Puffer	Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Quelle	E.coli
Reinheit	Greater than 85% as determined by SDS-PAGE.
Formulierung	Lyophilized powder

Sequenz	MTGPHSNIVH VIGAGLAGSE AAWQVAKAGV PVMLHEMRPD RMTEAHRTDG LAELVCSNSF RSDDAANNAV GLLHAEMRRL GSLIMRAADA NQVPAGGALA VDRDGFSAAV TKALNDHPLI EIARGEIAGL PPADWGNVIV ATGPLTSAPL ADAIRELTDE NALAFFDAIA PIVHRESIDM SVAWFQSRDYD KVGPGGNGAD YINCPMTKEQ YEGFVAALIA GEKTEFKEWE TN
Anwendungsbeschreibung	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.