

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Influenza A virus Nucleoprotein antibody [C43] (HRP), IgG2a, Mouse, Monoclonal GTX00843

|                          |                                                                                |
|--------------------------|--------------------------------------------------------------------------------|
| Artikelname              | Influenza A virus Nucleoprotein antibody [C43] (HRP), IgG2a, Mouse, Monoclonal |
| Artikelnummer            | GTX00843                                                                       |
| Hersteller Artikelnummer | GTX00843                                                                       |
| Alternativnummer         | GTX00843-50                                                                    |
| Hersteller               | GeneTex                                                                        |
| Wirt                     | Mouse                                                                          |
| Kategorie                | Antikörper                                                                     |
| Applikation              | ELISA, ICC, WB                                                                 |
| Spezies Reaktivität      | Virus                                                                          |
| Immunogen                | Human Influenza A Virus H3N2                                                   |
| Konjugation              | HRP                                                                            |
| Klonalität               | Monoclonal                                                                     |
| Konzentration            | 1 mg/ml (Please refer to the vial label for the specific concentration.)       |
| Klon-Bezeichnung         | [C43]                                                                          |
| Isotyp                   | IgG2a                                                                          |

|                        |                                                                                                                                                                                                                                                                 |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Sensitivitaet          | It reacts with NP of all influenza A viruses so far tested, including seasonal H2N2, H3N2(A/Sydney/5/1997), and H5N1(A/crow/Kyoto53/2004), H5N1 (A/duck/Egypt/D2br10/07), H5N1(A/duck/HK/342/78), H5N2(A/crow/Kyoto/53/04), H9N1, H9N2 (A/Turkey/Wisconsin/1/66 |
| Puffer                 | Filter-sterilized PBS, 50% Glycerol, no preservative.                                                                                                                                                                                                           |
| Reinheit               | Purified IgG                                                                                                                                                                                                                                                    |
| Formulierung           | Liquid                                                                                                                                                                                                                                                          |
| Anwendungsbeschreibung | WB: ~1:1000. ICC/IF: ~1:200. *Optimal dilutions/concentrations should be determined by the researcher. Not tested in other applications.                                                                                                                        |