

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Trident RIPA Lysis Buffer GTX400005

Artikelname	Trident RIPA Lysis Buffer
Artikelnummer	GTX400005
Hersteller Artikelnummer	GTX400005
Alternativnummer	GTX400005-100
Hersteller	GeneTex
Kategorie	Sonstiges
Applikation	WB
Produktbeschreibung	Radio-Immunoprecipitation assay lysis buffer The ingredients of RIPA lysis buffer are as follows: 25 mM Tris pH 7.6 150 mM NaCl 1% Triton X-100 1% sodium deoxycholate 1% SDS...
Formulierung	Liquid
Anwendungsbeschreibung	This product contains 1% SDS, which is higher than what is in some other commercial RIPA lysis buffers, and should facilitate a more thorough extraction of nuclear protein. The sample may be viscous following extraction due to unsheared DNA. To decrease this viscosity, we recommend increasing the volume of GTX400005 and trying either sonication or drawing the sample back-and-forth gently (to avoid bubbles) ten times through an 18G blunt needle attached to a 1ml syringe.