

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Neutrophil Alkaline Phosphatase (NAP) ELISA KBS-KBH1950

|                          |                                                                                                                                                                                                                                                                |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Neutrophil Alkaline Phosphatase (NAP) ELISA                                                                                                                                                                                                      |
| Artikelnummer            | KBS-KBH1950                                                                                                                                                                                                                                                    |
| Hersteller Artikelnummer | KBH1950                                                                                                                                                                                                                                                        |
| Alternativnummer         | KBS-KBH1950-1                                                                                                                                                                                                                                                  |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                            |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                    |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                          |
| Produktbeschreibung      | The GENLISA Human Neutrophil Alkaline Phosphatase (NAP) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validated against... |
| Detektionsbereich        | 1.563 - 100 ng/ml                                                                                                                                                                                                                                              |
| Sensitivität             | 0.75 ng/ml                                                                                                                                                                                                                                                     |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                             |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                     |