

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### GENLISA Human Trophoblast Antigen-Immunoglobulin G (TA-IgG) ELISA KBS-KBH3568

|                          |                                                                                                                                                                                                                                                               |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Artikelname              | GENLISA Human Trophoblast Antigen-Immunoglobulin G (TA-IgG) ELISA                                                                                                                                                                                             |
| Artikelnummer            | KBS-KBH3568                                                                                                                                                                                                                                                   |
| Hersteller Artikelnummer | KBH3568                                                                                                                                                                                                                                                       |
| Alternativnummer         | KBS-KBH3568-1                                                                                                                                                                                                                                                 |
| Hersteller               | Krishgen Biosystems                                                                                                                                                                                                                                           |
| Kategorie                | Kits/Assays                                                                                                                                                                                                                                                   |
| Spezies Reaktivität      | Human                                                                                                                                                                                                                                                         |
| Produktbeschreibung      | The GENLISA Human Trophoblast Antigen-Immunoglobulin G (TA-IGG) ELISA includes features like: Ready to use protocol with break-apart wells for ease of use Standardisation and High Reproducibility Lot to Lot Consistency Accuracy and Precision Validate... |
| Detektionsbereich        | Qualitative                                                                                                                                                                                                                                                   |
| Sensitivitaet            | ---                                                                                                                                                                                                                                                           |
| Proben                   | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. For tear samples, please connect with us to discuss further.                                                                                            |
| Anwendungsbeschreibung   | Colorimetric, Absorbance measured at 450nm                                                                                                                                                                                                                    |