

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Protein Molecular Weight Marker, Unconjugated Preis auf Anfrage BYT-ORB1784560

Artikelname	Protein Molecular Weight Marker, Unconjugated Preis auf Anfrage
Artikelnummer	BYT-ORB1784560
Hersteller Artikelnummer	orb1784560
Alternativnummer	BYT-ORB1784560-5
Hersteller	Biorbyt
Kategorie	Proteine/Peptide
Applikation	SDS-PAGE
Konjugation	Unconjugated
Produktbeschreibung	Protein Molecular Weight Marker...
Konzentration	1X
Puffer	Preservative: 0.01% (w/v) Sodium Azide. Stabilizer: None. 2.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) beta-Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue
Reinheit	Protein Molecular Weight Marker is supplied in 1X loading buffer and is ready-to-use. It contains: [62.5 mM Tris-HCl, pH 6.8, 2% (w/v) SDS, 250 mM NaCl, 1mM EDTA, 4% (v/v) beta-Mercaptoethanol, 35% (v/v) Glycerol and 0.05% (w/v) bromophenol blue.] This product are unstained markers optimized for use with Coomassie blue staining, but the loading volume can be adjusted for compatibility with other staining methods (i.e. silver staining). The intensity of the 67 kDa bovine albumin band has been increased to serve as a reference.

Formulierung	Liquid (in 1x Loading Buffer)
Anwendungsbeschreibung	Application Notes: Protein Molecular Weight Marker has been tested and is designed for use in accurate estimation of protein molecular weight in SDS-polyacrylamide gel electrophoresis and western blotting. Thaw the cocktail at room temperature. Gently invert vial several times to ensure that the solution is homogeneous and that any precipitated material is re-dissolved. A loading volume of 5-10µl is appropriate for most gel formats