

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Tris-Glycine-Native Electrophoresis Buffer Powder ABB-BR00006

Artikelname	Tris-Glycine-Native Electrophoresis Buffer Powder
Artikelnummer	ABB-BR00006
Hersteller Artikelnummer	BR00006
Alternativnummer	ABB-BR00006-10X500ML, ABB-BR00006-5X500ML
Hersteller	ABclonal
Kategorie	Sonstiges
Produktbeschreibung	Tris-Glycine-Native Electrophoresis Buffer is compatible with Tris-based gels and is commonly used for separating large-molecular-weight proteins. It is designed for native polyacrylamide gel electrophoresis (Native-PAGE).This product is supplied in ...
Konzentration	Protein electrophoresis
Application Verdünnung	Product CodeProduct NameSpecificationStorage Temperature BR00006Tris-Glycine-Native Electrophoresis Buffer Powder5*500 mL / 10*500 mL Store at room temperature
Anwendungsbeschreibung	Cross-Reactivity: 1.Fast. 2.Clear protein bands. 3.Simple operation,no need to weigh or filter reagents.