

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Recombinant SARS-CoV-2 Spike S1+S2 ECD(S-ECD)(B.1.1.529/Omicron) Trimer Protein, Human ABB-RP02111

Artikelname	Recombinant SARS-CoV-2 Spike S1+S2 ECD(S-ECD)(B.1.1.529/Omicron) Trimer Protein, Human
Artikelnummer	ABB-RP02111
Hersteller Artikelnummer	RP02111
Alternativnummer	ABB-RP02111-100UG
Hersteller	ABclonal
Wirt	Human
Kategorie	Proteine/Peptide
Spezies Reaktivität	Virus
Immunogen	Val16-Pro1213
Produktbeschreibung	The spike protein (S) of coronavirus (CoV) attaches the virus to its cellular receptor, angiotensin-converting enzyme 2 (ACE2). A defined receptor-binding domain (RBD) on S mediates this interaction. The S protein plays key parts in the induction of n...
Konzentration	< 1 EU/µg of the protein by LAL method.
Molekulargewicht	136.70 kDa
NCBI	43740568
UniProt	O60594
Reinheit	95 % as determined by SDS-PAGE.

Formulierung	Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Target-Kategorie	SARS-CoV-2 Spike S1+S2 ECD(S-ECD)(B.1.1.529/Omicron) Trimer
Application Verdünnung	Lyophilized from a 0.22 µm filtered solution of PBS, pH 7.4.Contact us for customized product form or formulation.
Anwendungsbeschreibung	Cross-Reactivity: Centrifuge the tube before opening. Reconstitute to a concentration of 0.1-0.5 mg/mL in sterile distilled water. Avoid vortex or vigorously pipetting the protein. For long term storage,it is recommended to add a carrier protein or stablizer (e.g. 0.1% BSA,5% HSA,10% FBS or 5% Trehalose),and aliquot the reconstituted protein solution to minimize free-thaw cycles., ResearchArea: SARS-CoV-2 antigens