

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

Product Datasheet

Rat NT-ProANP(N-Terminal Pro-Atrial Natriuretic Peptide) Microsample ELISA Kit EBT-ELK2400MS

Artikelname	Rat NT-ProANP(N-Terminal Pro-Atrial Natriuretic Peptide) Microsample ELISA Kit
Artikelnummer	EBT-ELK2400MS
Hersteller Artikelnummer	ELK2400MS
Alternativnummer	EBT-ELK2400MS-96, EBT-ELK2400MS-48, EBT-ELK2400MS-96X5
Hersteller	ELK Biotechnology
Kategorie	Kits/Assays
Spezies Reaktivität	Rat
Konzentration	40 ng/mL
Detektionsbereich	0.63-40 ng/mL
Sensitivitaet	0.246 ng/mL
Proben	Serum, plasma, tissue homogenates and other biological fluids.

Anwendungsbeschreibung

Assay Type: Sandwich. Assay length: 3.5h. Research Area: Endocrinology, Cardiovascular biology,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rat NT-ProANP. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rat NT-ProANP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rat NT-ProANP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Rat NT-ProANP in the samples is then determined by comparing the OD of the samples to the standard curve