

Bitte beachten Sie: Dieses Dokument wurde automatisch erstellt und ist kein Ersatz für das Originaldokument des Herstellers.

## Product Datasheet

### Mouse CCR2(Chemokine C-C-Motif Receptor 2) ELISA Kit EBT-ELK5338

Artikelname	Mouse CCR2(Chemokine C-C-Motif Receptor 2) ELISA Kit
Artikelnummer	EBT-ELK5338
Hersteller Artikelnummer	ELK5338
Alternativnummer	EBT-ELK5338-96, EBT-ELK5338-48, EBT-ELK5338-96X5
Hersteller	ELK Biotechnology
Kategorie	Kits/Assays
Spezies Reaktivität	Mouse
Konzentration	5000 pg/mL
Detektionsbereich	78.13-5000 pg/mL
Sensitivitaet	33 pg/mL
UniProt	<a href="#">P51683</a>
Proben	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Anwendungsbeschreibung	<p>Assay Type: Sandwich. Assay length: 3.5h. Research Area: CD &amp; Adhesion molecule, Tumor immunity, Infection immunity,. Test principle: The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Mouse CCR2. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Mouse CCR2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Mouse CCR2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm 10nm. The concentration of Mouse CCR2 in the samples is then determined by comparing the OD of the samples to the standard curve</p>
------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------